



TspDTI

TspDT I

Restriction Endonuclease

Recognition Sequence:



Cat. No.	Size
E2502-01	50 units
E2502-02	250 units

Reaction Temperature: 70°C

Inactivation Temperature (20 min): --

Prototype: TspDTI

Source: *Thermus species* DT

Note 1: Purified from *E.coli* strain that carries the cloned tspDTRI gene from *Thermus* sp. DT.

Package Contents:

- TspDTI
- 10x Reaction Buffer TspDTI

Storage Conditions: Store at -20°C

Recommended Buffer: TspDTI
(or compatible third party buffers)

DNA Methylation:

No inhibition: dam, dcm, EcoKI, CpG

Standard Reaction Protocol:

Mix the following reaction components:

- 1-2 µg pure DNA or 10 µl PCR product (= -0.1-2 µg DNA)
- 3 µl 10x Buffer TspDTI

1-2 U TspDTI (use 1 U / µg DNA, < 10 % React. Volume!)

Tips: Add enzyme as last component. Mix components well before adding enzyme. After enzyme addition, mix gently by pipetting. Do not vortex.

@ 30 µl H₂O, DNA and DNase free

Incubate for 3 h at 70°C

Stop reaction by alternatively

- (a) Addition of 1.2 µl EDTA pH 8.0 (0.5 M), final 20 mM or
- (b) Heat Inactivation

20 min at 89°C (not recommended) or

- (c) Spin Column DNA Purification
(e.g. EURx PCR/DNA CleanUp Kit, Cat.No. E3520) or
- (d) Gel Electrophoresis and Single Band Excision
(e.g. EURx AgaroseOut DNA Kit, Cat.No. E3540) or
- (e) Phenol-Chloroform Extraction or Ethanol Precipitation.

Note 1: It is required to purify DNA before digestion. We recommend PCR/DNA Clean-Up Purification (E3520) or Agarose Out DNA Purification Kit (E3540).

Note 2: It is not recommended to use more than 2 units per 30 µl reaction. It is strongly suggested to perform digestion for over 1 hr.

Note 3: To avoid DNA shift during electrophoresis caused by strong protein-DNA interaction, it is recommended to terminate reaction by addition of reaction stop solution (containing denaturing reagent, i.e. 0,2% [w/v] SDS) followed by 20 minutes heat inactivation at 89°C.

Unit Definition:

One unit is the amount of enzyme required to digest 1 µg of pUC19 DNA to obtain stable digestion pattern in 1 hr in a total reaction volume of 30 µl. Enzyme activity was determined in the recommended reaction buffer.

Reaction Buffer:

1 x TspDTI Buffer: 10 mM Tris-HCl (pH 8.5 at 25°C), 10 mM MgCl₂, 1 mM dithiothreitol + enhancers (1).

Avoid multiple cycles of freezing/thawing of the stock reaction buffer /no more than 3 times/. Thawing should be performed at temperatures not exceeding 10°C. Recommended procedure is to divide the provided reaction buffer into smaller portions and preserve them at -70°C for long-term. Temperature of -20°C should be used only for short-term storage.

Storage Buffer:

20 mM Tris-HCl, (pH 8.3 at 25°C), 25 mM (NH₄)₂SO₄, 25 mM KCl, 0,5 mM EDTA, 0,5 mM dithiothreitol, 0,02 % [v/v] Tergitol™ TMN, 0,02 % Tween™20, 0,02 % Igepal, 50 % [v/v] glycerol.

Quality Control:

All preparations are assayed for contaminating endonuclease, 3'-exonuclease, 5'-exonuclease/5'-phosphatase, as well as nonspecific single- and double-stranded DNase activities.

References:

1. Skowron, P.M., Majewski, J., Żylicz-Stachula, A., Rutkowska, S.M., Jaworowska, I. and Harasimowicz-Słowińska, R.I. (2003). Nucleic Acids Research 31, 14 e 74.