

Perfect 25 Protein Ladder

Perfect 25 Protein Ladder

Protein ladder, composed of precisely spaced recombinant bands, for sizing polypeptides ranging from 25 kDa to 150 kDa.

Cat. No.**Size****Description:**

E3200-01

100 loadings

→ The only protein ladder available in precise increments with 7 bands ranging from 25 to 150 kDa.

E3200-02

500 loadings

→ Comes in a ready-to-use format.

Storage Conditions:

Store at -20°C

Loading:

The recommended amount of size marker to load on a gel is 5-10 µl per lane depending on a gel type and size of well.

Concentration:

Concentration is approximately 1 mg/ml and is supplied in gel loading buffer.

